

Sample Name :C-Type Natriuretic Peptide (CNP) (1-22), human
Cat No. :RP11110
Lot No. :E9851621I
Time Processed : 10:09:35
Month-Day-Year Processed : 11/07/2021

Pump A : 0.065% trifluoroacetic in 100% water (v/v)
Pump B : 0.05% trifluoroacetic in 100% acetonitrile (v/v)
Total Flow:1 ml/min
Wavelength:220 nm

Time	Module	Command	Value
0.01	Pumps	B.Conc	5
25.00	Pumps	B.Conc	65
25.01	Pumps	B.Conc	95
27.00	Pumps	B.Conc	95
27.01	Pumps	B.Conc	5
35.00	Pumps	B.Conc	5
35.01	Controller	Stop	


<<Column Performance>>

<Detector A>

Column : Inertsil ODS-3 4.6 x 250 mm

Equipment: GK11010017

mV



1 Detector A Channel 1 / 220nm

Peak Table

Detector A Channel 1 220nm

Peak#	Ret. Time	Area	Height	Area %
1	12.613	12400	556	0.275
2	16.067	8315	1132	0.184
3	17.232	30067	3518	0.666
4	17.721	42618	2411	0.945
5	18.090	4355453	571989	96.534
6	18.563	43878	4619	0.973
7	18.797	17269	1547	0.383
8	20.055	1818	218	0.040
Total		4511818	585991	100.000