


Mass Spectrum



Sample Information

User :	Gary	Equipment	: GK11010007
Month-Day Processed :	02/28/20	Interface Bias	: +4.5kV
Time Processed :	02:50:55 PM	Drying Gas Flow	: 5 L/min
Injection Volume :	0.5	T.Flow	: 0.2ml/min
Sample Name :	Apidaecin IB	B.conc	: 50%H2O/50%MeOH
Cat No. :	RP20409		
Lot No. :	E2862620B		
Theoretical MW :	2108.40		
Observed MW :	2108.0		