


Mass Spectrum



Sample Information

Month-Day Processed : 11/26/21
Time Processed : 10:41:03 AM
Injection Volume : 0.4
Sample Name : TAT peptide
Cat No. : RP20256
Lot No. : E3466921K
Theoretical MW : 1621.91
Observed MW : 1622.0

Interface : ESI
Nebulizing Gas Flow : 1.5L/min
CDL Temp : 250
Block Temp : 200

Equipment : GK11010007
Interface Bias : +3.5 kV
Drying Gas Flow : 5 L/min
T.Flow : 0.2 ml/min
B.conc : 50%H₂O/50%MeOH