


Mass Spectrum



Sample Information

Month-Day Processed : 08/23/20
 Time Processed : 0:01:55
 Injection Volume : 0.3
 Sample Name : Tetanus toxin (830-843)
 Cat. No. : RP19980
 Lot. No. : E6884020H
 Theoretical MW : 1611.84
 Observed MW : 1611.6

Interface : ESI

Nebulizing Gas Flow : 1.5L/min
 CDL Temp : 250
 Block Temp : 200

Equipment

Interface Bias : +4.5 kV
 Drying Gas Flow : 5 L/min
 T.Flow : 0.2 ml/min
 B.conc : 50% H₂O/50% MeOH

Equipment : GK11010007