


Mass Spectrum



Sample Information

Month-Day Processed : 12/26/20
Time Processed : 10:15:42 PM
Injection Volume : 0.4
Sample Name : Myomodulin
Cat No. : RP10606
Lot No. : E4291520L
Theoretical MW : 846.20
Observed MW : 845.8

Interface : ESI
Nebulizing Gas Flow : 1.5L/min
CDL Temp : 250
Block Temp : 200

Equipment : GK11010007
Interface Bias : +3.5 kV
Drying Gas Flow : 5 L/min
T.Flow : 0.2 ml/min
B.conc : 50%H₂O/50%MeOH