


Mass Spectrum



Sample Information

Acquired by : Gary
Month-Day Processed : 10/06/18
Time Processed : 14:47:49
Injection Volume : 0.2
Sample Name : HIF-1 alpha (556-574)
Cat. No. : RP20479
Lot. No. : PE2891809
Theoretical MW : 2254.60
Observed MW : 2255.0

Interface : ESI
Nebulizing Gas Flow : 1.5L/min
CDL Temp : 250°C
Block Temp : 200°C

Equipment : GK11010007
Interface Bias : +4.5 kV
Drying Gas Flow : 5 L/min
T.Flow : 0.2 ml/min
B.conc : 50%H₂O/50%MeOH