


Mass Spectrum



Sample Information

Month-Day Processed : 01/27/21
Time Processed : 4:44:22 AM
Injection Volume : 0.6
Sample Name : HIF-1 alpha (556-574)
Cat No. : RP20479
Lot No. : E0645321A
Theoretical MW : 2254.60
Observed MW : 2254.0

Interface : ESI
Nebulizing Gas Flow : 1.5L/min
CDL Temp : 250
Block Temp : 200

Equipment : GK11010007
Interface Bias : +3.5 kV
Drying Gas Flow : 5 L/min
T.Flow : 0.2 ml/min
B.conc : 50%H₂O/50%MeOH