


Mass Spectrum



Sample Information

Month-Day Processed : 11/28/20
Time Processed : 10:33:50 PM
Injection Volume : 0.3
Sample Name : Brain Natriuretic Peptide (BNP) (1-32), human
Cat No. : RP111119
Lot No. : E2728420K
Theoretical MW : 3464.10
Observed MW : 3463.0

Interface : ESI
Nebulizing Gas Flow : 1.5L/min
CDL Temp : 250
Block Temp : 200

Equipment : GK11010007
Interface Bias : +3.5 kV
Drying Gas Flow : 5 L/min
T.Flow : 0.2 ml/min
B.conc : 50%H2O/50%MeOH