


Mass Spectrum



Sample Information

User : Gary
Month-Day Processed : 03/18/19
Time Processed : 07:53:38 PM
Injection Volume : 0.2
Sample Name : YRGDS Fibronectin Fragment
Cat No. : RP20476
Lot No. : PE9291903
Theoretical MW : 596.60
Observed MW : 596.6

Interface : ESI
Nebulizing Gas Flow : 1.5L/min
CDL Temp : 250
Block Temp : 200

Equipment : GK11010007
Interface Bias : +4.5kV
Drying Gas Flow : 5 L/min
T.Flow : 0.2ml/min
B.conc : 50%H₂O/50%MeOH