


Mass Spectrum



Sample Information

User :	Gary	Equipment	: GK11010007
Month-Day Processed :	07/13/20	Interface Bias	: +4.5kV
Time Processed :	03:23:29 PM	Drying Gas Flow	: 5 L/min
Injection Volume :	0.3	T.Flow	: 0.2ml/min
Sample Name :	Fibronectin Active Fragment Control	B.conc	: 50% H2O/50% MeOH
Cat No. :	RP20347		
Lot No. :	E2308720G		
Theoretical MW :	504.50		
Observed MW :	504.4		