


# Mass Spectrum



Sample Information  
Month-Day Processed : 01/07/21  
Time Processed : 9:42:20 AM  
Injection Volume : 0.3  
Sample Name : LL37 (Human)  
Cat No. : RP13323  
Lot No. : E8052221A  
Theoretical MW : 4493.28  
Observed MW : 4492.4

Interface : ESI  
Nebulizing Gas Flow : 1.5L/min  
CDL Temp : 250  
Block Temp : 200

Equipment : GK11010007  
Interface Bias : +3.5 kV  
Drying Gas Flow : 5 L/min  
T.Flow : 0.2 ml/min  
B.conc : 50%H2O/50%MeOH