


Mass Spectrum



Sample Information		Interface		Equipment	
Month-Day Processed :	10/01/21	Nebulizing Gas Flow : 1.5L/min		: GK11010007	
Time Processed :	7:29:45	CDL Temp : 250		: +3.5 kV	
Injection Volume :	0.2	Block Temp : 200		: 5 L/min	
Sample Name :	DYKDDDDK Peptide			: 0.2 ml/min	
Cat No. :	RP10586			: 50% H ₂ O/50% MeOH	
Lot No. :	E0319621I				
Theoretical MW :	1012.97				
Observed MW :	1013.0				