

Sample Name :Glucagon-Like Peptide (GLP) II, human
 Cat No. :RP10769
 Lot No. :E3439221K
 Time Processed :9:16:11
 Month-Day-Year Processed :11/27/2021


Pump A : 0.065% trifluoroacetic in 100% water (v/v)
 Pump B : 0.05% trifluoroacetic in 100% acetonitrile (v/v)
 Total Flow:1 ml/min
 Wavelength:220 nm

Time	Module	Command	Value
0.01	Pumps	B.Conc	5
25.00	Pumps	B.Conc	65
25.01	Pumps	B.Conc	95
27.00	Pumps	B.Conc	95
27.01	Pumps	B.Conc	5
35.00	Pumps	B.Conc	5
35.01	Controller	Stop	

<<Column Performance>>

<Detector A>

Column : Inertsil ODS-3 4.6 x 250 mm



1 Detector A Channel 1 / 220nm

Peak Table

Detector A Channel 1 220nm

Peak#	Ret. Time	Area	Height	Area %
1	16.026	42181	8153	0.255
2	16.494	128315	15004	0.774
3	16.600	48063	9945	0.290
4	16.779	61570	8405	0.371
5	17.016	15905616	2308500	95.970
6	17.366	336394	39721	2.030
7	17.600	51363	8849	0.310
Total		16573503	2398577	100.000