


Mass Spectrum



Sample Information
Month-Day Processed : 03/18/21
Time Processed : 9:23:26 AM
Injection Volume : 0.4
Sample Name : LL-37, acetylated, amidated
Cat No. : RP20506
Lot No. : E1292321C
Theoretical MW : 4534.40
Observed MW : 4534.8

Interface : ESI
Nebulizing Gas Flow : 1.5L/min
CDL Temp : 250
Block Temp : 200

Equipment : GK11010007
Interface Bias : +3.5 kV
Drying Gas Flow : 5 L/min
T.Flow : 0.2 ml/min
B.conc : 50%H2O/50%MeOH