


Mass Spectrum



Sample Information

Acquired by : Gary
Month-Day Processed : 06/05/19
Time Processed : 08:06:22 PM
Injection Volume : 0.3
Sample Name : Brain Natriuretic Peptide (BNP) (1-32), human
Cat. No. : RP11119
Lot. No. : PE1781905
Theoretical MW : 3464.10
Observed MW : 3463.8

Interface

Interface : ESI
Nebulizing Gas Flow : 1.5L/min
CDL Temp : 250
Block Temp : 200

Equipment

Equipment : GK11010007
Interface Bias : +4.5 kV
Drying Gas Flow : 5 L/min
T.Flow : 0.2 ml/min
B.conc : 50%H2O/50%MeOH