

Sample Name :Glucagon-Like Peptide (GLP) II, human
 Cat No. :RP10769
 Lot No. :E1884721L
 Time Processed :22:04:57
 Month-Day-Year Processed :12/08/2021


Pump A : 0.065% trifluoroacetic in 100% water (v/v)
 Pump B : 0.05% trifluoroacetic in 100% acetonitrile (v/v)
 Total Flow:1 ml/min
 Wavelength:220 nm

Time	Module	Command	Value
0.01	Pumps	Pump A B.Conc	5
25.00	Pumps	Pump A B.Conc	65
25.01	Pumps	Pump A B.Conc	95
27.00	Pumps	Pump A B.Conc	95
27.01	Pumps	Pump A B.Conc	5
35.00	Pumps	Pump A B.Conc	5
35.01	Controller	Stop	

<<Column Performance>>

<Detector A>

Column : Inertsil ODS-3 4.6 x 250 mm



1 Detector A Channel 1 / 220nm

Peak Table

Detector A Channel 1 220nm

Peak#	Ret. Time	Area	Height	Area %
1	13.236	21160	3022	0.179
2	13.442	1815	324	0.015
3	13.542	1465	343	0.012
4	13.802	96214	10593	0.814
5	14.484	1210	271	0.010
6	14.670	2360	290	0.020
7	14.785	1500	312	0.013
8	14.986	3507	547	0.030
9	15.351	6240	610	0.053
10	15.591	1691	288	0.014
11	15.929	27634	2936	0.234
12	16.458	8383	1178	0.071
13	16.735	2181	439	0.018
14	16.927	30621	4148	0.259
15	17.619	56780	3307	0.481
16	17.908	11357937	877293	96.124
17	18.440	188515	18146	1.595
18	19.280	1066	229	0.009
19	19.519	5595	715	0.047
Total		11815874	924990	100.000