


Mass Spectrum



Sample Information

Acquired by : Gary

Month-Day Processed : 10/19/17

Time Processed : 17:10:34

Injection Volume : 0.2

Sample Name : Human Papillomavirus (HPV) E7 protein (49 - 57)

EcVVP q0².....TR4246:

NqVVP q0².....RC9983932

Theoretical MW : 1120.27

Observed MW : 1120.30

Interface : ESI

Nebulizing Gas Flow : 1.5L/min

CDL Temp : 250°C

Block Temp : 200°C

Interface Bias : +4.5 kV

Drying Gas Flow : 5 L/min

T.Flow : 0.2 ml/min

B.conc : 50% H₂O/50% MeOH