


Mass Spectrum



Sample Information

User : Gary
 Month-Day Processed : 03/02/19
 Time Processed : 01:43:17 PM
 Injection Volume : 0.2
 Sample Name : BNP (1-32), human
 Cat No. : RP11119
 Lot No. : PE7971902
 Theoretical MW : 3464.10
 Observed MW : 3463.8

Interface : ESI
 Nebulizing Gas Flow : 1.5L/min
 CDL Temp : 250
 Block Temp : 200

Equipment : GK11010007
 Interface Bias : +4.5kV
 Drying Gas Flow : 5 L/min
 T.Flow : 0.2ml/min
 B.conc : 50%H2O/50%MeOH