


Mass Spectrum



Sample Information

Month-Day Processed : 10/23/20
Time Processed : 21:27:26
Injection Volume : 0.3
Sample Name : Proctolin
Cat. No. : RP10475
Lot. No. : E2527620J
Theoretical MW : 648.75
Observed MW : 648.3

Interface : ESI
Nebulizing Gas Flow : 1.5L/min
CDL Temp : 250
Block Temp : 200

Equipment : GK11010007
Interface Bias : +4.5 kV
Drying Gas Flow : 5 L/min
T.Flow : 0.2 ml/min
B.conc : 50% H2O/50% MeOH