


Mass Spectrum



Sample Information	Interface	Interface Bias
Acquired by : Gary	Nebulizing Gas Flow : 1.5L/min	Drying Gas Flow : +4.5 kV
Month-Day Processed : 05/05/17	CDL Temp : 250°C	T.Flow : 5 L/min
Time Processed : 15:17:29	Block Temp : 200°C	B.conc : 0.2 ml/min
Injection Volume : 0.2		
Sample Name : Human Papillomavirus (HPV) E7 protein (49 - 57)		
Cat. No. : RP20249		
Lot. No. : PE6801704		
Theoretical MW : 1120.27		
Observed MW : 1120.40		