


Mass Spectrum



Sample Information

User : Gary
 Month-Day Processed : 03/02/19
 Time Processed : 08:14:40 PM
 Injection Volume : 0.2
 Sample Name : LL-37, scrambled
 Cat No. : RP20332
 Lot No. : PE8951902
 Theoretical MW : 4493.30
 Observed MW : 4492.6

Interface : ESI
 Nebulizing Gas Flow : 1.5L/min
 CDL Temp : 250
 Block Temp : 200

Equipment : GK11010007
 Interface Bias : +4.5kV
 Drying Gas Flow : 5 L/min
 T.Flow : 0.2ml/min
 B.conc : 50%H2O/50%MeOH