


Mass Spectrum



Sample Information

User : Gary
Month-Day Processed : 04/04/20
Time Processed : 10:09:12 AM
Injection Volume : 0.3
Sample Name : Somatostatin
Cat No. : RP10226
Lot No. : E6035620C
Theoretical MW : 1637.88
Observed MW : 1637.4

Interface : ESI
Nebulizing Gas Flow : 1.5L/min
CDL Temp : 250
Block Temp : 200

Equipment : GK11010007
Interface Bias : +4.5kV
Drying Gas Flow : 5 L/min
T.Flow : 0.2ml/min
B.conc : 50%H2O/50%MeOH