


Mass Spectrum



Sample Information
Month-Day Processed : 12/25/20
Time Processed : 11:49:07 AM
Injection Volume : 0.4
Sample Name : MDYKDHGDYKDHDDYKDDDDK peptide
Cat No: RP21087
Lot No: E4936320L
Theoretical MW : 2861.87
Observed MW : 2860.8

Interface : ESI
Nebulizing Gas Flow : 1.5L/min
CDL Temp : 250
Block Temp : 200

Equipment : GK11010007
Interface Bias : +4.5 kV
Drying Gas Flow : 5 L/min
T.Flow : 0.2 ml/min
B.conc : 50% H₂O/50% MeOH