


# Mass Spectrum



Sample Information  
Month-Day Processed : 09/11/20  
Time Processed : 03:08:45 PM  
Injection Volume : 0.4  
Sample Name : Gastrin-1, human  
Cat. No. : RP12740  
Lot. No. : E18497201  
Theoretical MW : 2098.21  
Observed MW : 2097.6

Interface : ESI  
Nebulizing Gas Flow : 1.5L/min  
CDL Temp : 250  
Block Temp : 200

Equipment : GK11010007  
Interface Bias : +4.5 kV  
Drying Gas Flow : 5 L/min  
T.Flow : 0.2 ml/min  
B.conc : 50%H2O/50%MeOH