


Mass Spectrum



Sample Information

Month-Day Processed : 12/22/20
Time Processed : 05:16:55 PM
Injection Volume : 0.3
Sample Name : Calcitonin Gene Related Peptide (CGRP), rat
Cat No : RP11095
Lot No : E7150020L
Theoretical MW : 3806.25
Observed MW : 3805.2

Interface : ESI
Nebulizing Gas Flow : 1.5L/min
CDL Temp : 250
Block Temp : 200

Equipment : GK11010007
Interface Bias : +4.5 kV
Drying Gas Flow : 5 L/min
T.Flow : 0.2 ml/min
B.conc : 50% H₂O/50% MeOH