


Mass Spectrum



Sample Information

Month-Day Processed : 04/21/21
Time Processed : 22:49:50
Injection Volume : 0.4
Sample Name : Influenza Hemagglutinin (307-319)
Cat No. : RP20087
Lot No. : E4296921D
Theoretical MW : 1503.79
Observed MW : 1503.6

Interface : ESI
Nebulizing Gas Flow : 1.5L/min
CDL Temp : 250
Block Temp : 200

Equipment : GK11010007
Interface Bias : +3.5 kV
Drying Gas Flow : 5 L/min
T.Flow : 0.2 ml/min
B.conc : 50%H2O/50%MeOH