


Mass Spectrum



Sample Information

Acquired by: Gary
Month-Day Processed: 07/17/17
Time Processed: 11:17:06
Injection Volume: 0.2
Sample Name: Tetanus toxin (830-843)
Sample ID: U7921CG060-1
Theoretical MW: 1611.84
Observed MW: 1611.5

Interface: ESI
Nebulizing Gas Flow: 1.5L/min
CDL Temp: 250°C
Block Temp: 200°C

Interface Bias
Drying Gas Flow
T.Flow
B.conc

: +4.5 kV
: 5 L/min
: 0.2 ml/min
: 50% H2O/50% MeOH