


Mass Spectrum



Sample Information

User : Gary
Month-Day Processed : 11/16/19
Time Processed : 9:00:12
Injection Volume : 0.3
Sample Name : Enfuvirtide (T-20)
Cat No. : RP17307
Lot No. : E6895119K
Theoretical MW : 4491.88
Observed MW : 4491.6

Interface

Interface : ESI
Nebulizing Gas Flow : 1.5L/min
CDL Temp : 250
Block Temp : 200

Equipment

Equipment : GK11010007
Interface Bias : +4.5kV
Drying Gas Flow : 5 L/min
T.Flow : 0.2ml/min
B.conc : 50%H2O/50%MeOH