


Mass Spectrum



Sample Information

User :	Gary	Equipment	: GK11010007
Month-Day Processed :	01/20/20	Interface Bias	: +4.5kV
Time Processed :	08:33:38 PM	Drying Gas Flow	: 5 L/min
Injection Volume :	0.4	T.Flow	: 0.2ml/min
Sample Name :	Angiotensin II, human	B.conc	: 50%H2O/50%MeOH
Cat No. :	RP11257		
Lot No. :	E6861520A		
Theoretical MW :	1046.18		
Observed MW :	1046.0		