


Mass Spectrum



Sample Information

User : Gary
 Month-Day Processed : 09/18/19
 Time Processed : 08:42:15 PM
 Injection Volume : 0.3
 Sample Name : Casein Kinase Substrates 3
 Cat No. : RP20214
 Lot No. : E4854319I
 Theoretical MW : 2131.27
 Observed MW : 2130.9

Interface : ESI
 Nebulizing Gas Flow : 1.5L/min
 CDL Temp : 250
 Block Temp : 200

Equipment : GK11010007
 Interface Bias : +4.5kV
 Drying Gas Flow : 5 L/min
 T.Flow : 0.2ml/min
 B.conc : 50%H2O/50%MeOH