


Mass Spectrum



Sample Information

Acquired by : Gary
Month-Day Processed : 04/26/18
Time Processed : 8:50:13
Injection Volume : 0.2
Sample Name : EGFRvIII peptide (PEPvIII)
Sample ID : U7986DD240-1
Theoretical MW : 1634.82
Observed MW : 1634.90

Interface : ESI
Nebulizing Gas Flow : 1.5L/min
CDL Temp : 250°C
Block Temp : 200°C

Interface Bias : +4.5 kV
Drying Gas Flow : 5 L/min
T.Flow : 0.2 ml/min
B.conc : 50% H2O/50% MeOH