


Mass Spectrum



Sample Information

User : Gary
Month-Day Processed : 06/30/17
Time Processed : 17:37:45
Injection Volume : 0.2
Sample Name : A-779
Cat No. : RP11249
Lot No. : PE6531706
Theoretical MW : 873.00
Observed MW : 872.90

Interface : ESI
Nebulizing Gas Flow : 1.5L/min
CDL Temp : 250°C
Block Temp : 200°C

Interface Bias : +4.5kV
Drying Gas Flow : 5 L/min
T.Flow : 0.2ml/min
B.conc : 50%H₂O/50%MeOH