


# Mass Spectrum



## Sample Information

Acquired by : Gary  
Month-Day Processed : 07/19/20  
Time Processed : 08:49:59 AM  
Injection Volume : 0.3  
Sample Name : Proctolin  
Sample ID : U809HFG150-1  
Theoretical MW : 648.75  
Observed MW : 648.6

Interface : ESI  
Nebulizing Gas Flow : 1.5L/min  
CDL Temp : 250  
Block Temp : 200

Equipment : GK11010007  
Interface Bias : +4.5 kV  
Drying Gas Flow : 5 L/min  
T.Flow : 0.2 ml/min  
B.conc : 50% H<sub>2</sub>O / 50% MeOH