


Mass Spectrum



Sample Information

User : Gary
 Month-Day Processed : 09/23/19
 Time Processed : 09:43:33 PM
 Injection Volume : 0.3
 Sample Name : Neuropeptide Y(1-36)
 Cat No. : RP20520
 Lot No. : E4942819I
 Theoretical MW : 4272.70
 Observed MW : 4272.0

Interface : ESI
 Nebulizing Gas Flow : 1.5L/min
 CDL Temp : 250
 Block Temp : 200

Equipment : GK11010007
 Interface Bias : +4.5kV
 Drying Gas Flow : 5 L/min
 T.Flow : 0.2ml/min
 B.conc : 50%H₂O/50%MeOH