


Mass Spectrum



Sample Information
Month-Day Processed : 07/30/21
Time Processed : 7:32:14
Injection Volume : 0.2
Sample Name : Exendin-4
Cat No. : RPI0874
Lot No. : E7807221G
Theoretical MW : 4186.66
Observed MW : 4186.4

Interface : ESI
Nebulizing Gas Flow : 1.5L/min
CDL Temp : 250
Block Temp : 200

Equipment : GK11010007
Interface Bias : +3.5 kV
Drying Gas Flow : 5 L/min
T.Flow : 0.2 ml/min
B.conc : 50%H2O/50%MeOH