


Mass Spectrum



Sample Information

User : Gary
Month-Day Processed : 02/10/20
Time Processed : 07:42:52 PM
Injection Volume : 0.3
Sample Name : LL37 (Human)
Cat No. : RP13323
Lot No. : E3585820B
Theoretical MW : 4493.28
Observed MW : 4493.4

Interface : ESI
Nebulizing Gas Flow : 1.5L/min
CDL Temp : 250
Block Temp : 200

Equipment : GK11010007
Interface Bias : +4.5kV
Drying Gas Flow : 5 L/min
T.Flow : 0.2ml/min
B.conc : 50%H2O/50%MeOH