


Mass Spectrum



Sample Information

Month-Day Processed : 12/02/20
Time Processed : 7:45:25 PM
Injection Volume : 0.3
Sample Name : C3a (70-77)
Cat No. : RP20511
Lot No. : E2686520K
Theoretical MW : 823.90
Observed MW : 823.5

Interface : ESI
Nebulizing Gas Flow : 1.5L/min
CDL Temp : 250
Block Temp : 200

Equipment : GK11010007
Interface Bias : +3.5 kV
Drying Gas Flow : 5 L/min
T.Flow : 0.2 ml/min
B.conc : 50%H2O/50%MeOH