


Mass Spectrum



Sample Information

Acquired by : Gary
Month-Day Processed : 08/01/18
Time Processed : 15:17:37
Injection Volume : 0.2
Sample Name : Somatostatin 28 (1-14)
Cat. No. : RP10230
Lot. No. : PE8741807
Theoretical MW : 1528.70
Observed MW : 1528.8

Interface : ESI
Nebulizing Gas Flow : 1.5L/min
CDL Temp : 250°C
Block Temp : 200°C

Interface Bias : +4.5 kV
Drying Gas Flow : 5 L/min
T.Flow : 0.2 ml/min
B.conc : 50% H₂O/50% MeOH