


Mass Spectrum



Sample Information
Month-Day Processed : 07/24/21
Time Processed : 11:01:08
Injection Volume : 0.2
Sample Name : Hexa-His
Cat. No. : RP11737
Lot. No. : E2270521G
Theoretical MW : 840.85
Observed MW : 840.8

Interface : ESI
Nebulizing Gas Flow : 1.5L/min
CDL Temp : 250
Block Temp : 200

Equipment : ZJ21010035
Interface Bias : +4.5 kV
Drying Gas Flow : 5 L/min
T.Flow : 0.2 ml/min
B.conc : 50% H2O/50% MeOH