

Sample Name :Glucagon-Like Peptide (GLP) II, human
 Cat No. :RP10769
 Lot No. :E4146621J
 Time Processed :2:57:52
 Month-Day-Year Processed :11/02/2021


Pump A : 0.065% trifluoroacetic in 100% water (v/v)
 Pump B : 0.05% trifluoroacetic in 100% acetonitrile (v/v)
 Total Flow:1 ml/min
 Wavelength:220 nm

Time	Module	Command	Value
0.01	Pumps	B.Conc	5
25.00	Pumps	B.Conc	65
25.01	Pumps	B.Conc	95
27.00	Pumps	B.Conc	95
27.01	Pumps	B.Conc	5
35.00	Pumps	B.Conc	5
35.01	Controller	Stop	

<<Column Performance>>

<Detector A>

Column : Inertsil ODS-3 4.6 x 250 mm



1 Detector A Channel 1 / 220nm

Peak Table

Detector A Channel 1 220nm

Peak#	Ret. Time	Area	Height	Area %
1	17.322	104648	17697	0.685
2	17.953	139210	11584	0.911
3	18.842	14566959	2315104	95.320
4	18.996	261206	92072	1.709
5	19.109	184902	33123	1.210
6	19.483	25188	4514	0.165
Total		15282112	2474095	100.000