


Mass Spectrum



Sample Information

User : Gary
Month-Day Processed : 04/30/19
Time Processed : 02:09:09 PM
Injection Volume : 0.2
Sample Name : Shepherdin (79-87)
Cat No. : RP20410
Lot No. : PE9871904
Theoretical MW : 949.10
Observed MW : 948.5

Interface : ESI
Nebulizing Gas Flow : 1.5L/min
CDL Temp : 250
Block Temp : 200

Equipment : GK11010007
Interface Bias : -3.5 kV
Drying Gas Flow : 5 L/min
T.Flow : 0.2 ml/min
B.conc : 50%H2O/50%MeOH