


Mass Spectrum



Sample Information

Acquired by : Gary
Month-Day Processed : 12/04/17
Time Processed : 15:37:09
Injection Volume : 0.2
Sample Name : Kassinin
Cat. No.: RP10461
Lot. No.: PE2051711
Theoretical MW : 1334.55
Observed MW : 1334.80

Interface : ESI
Nebulizing Gas Flow : 1.5L/min
CDL Temp : 250°C
Block Temp : 200°C

Interface Bias : +4.5 kV
Drying Gas Flow : 5 L/min
T.Flow : 0.2 ml/min
B.conc : 50% H₂O/50% MeOH