


Mass Spectrum



Sample Information

Month-Day Processed : 12/26/20
Time Processed : 9:39:18 PM
Injection Volume : 0.4
Sample Name : Astressin
Cat No. : RP10980
Lot No. : E4293720L
Theoretical MW : 3563.16
Observed MW : 3562.0

Interface : ESI
Nebulizing Gas Flow : 1.5L/min
CDL Temp : 250
Block Temp : 200

Equipment : GK11010007
Interface Bias : +3.5 kV
Drying Gas Flow : 5 L/min
T.Flow : 0.2 ml/min
B.conc : 50%H2O/50%MeOH