


Mass Spectrum



Sample Information

Month-Day Processed : 11/27/20
Time Processed : 10:07:01 PM
Injection Volume : 0.3
Sample Name : Angiotensin II, human
Cat No. : RP11257
Lot No. : E5929120K
Theoretical MW : 1046.18
Observed MW : 1045.8

Interface : ESI
Nebulizing Gas Flow : 1.5L/min
CDL Temp : 250
Block Temp : 200

Equipment : GK11010007
Interface Bias : +3.5 kV
Drying Gas Flow : 5 L/min
T.Flow : 0.2 ml/min
B.conc : 50% H₂O/50% MeOH