


# Mass Spectrum



Sample Information  
Month-Day/Processed : 11/26/20  
Time Processed : 04:45:07 PM  
Injection Volume : 0.4  
Sample Name : EGFRvIII peptide (PEPvIII)  
Cat No. : RP20356  
Lot No. : E5938520K  
Theoretical MW : 1634.82  
Observed MW : 1634.2

Interface : ESI  
Nebulizing Gas Flow : 1.5L/min  
CDL Temp : 250  
Block Temp : 200

Equipment : GK11010007  
Interface Bias : +3.5 kV  
Drying Gas Flow : 5 L/min  
T.Flow : 0.2 ml/min  
B.conc : 50%H2O/50%MeOH