


Mass Spectrum



Sample Information

Acquired by : Gary
Month-Day Processed : 03/28/18
Time Processed : 15:19:55
Injection Volume : 0.2
Sample Name : SLLK, Control Peptide for TSP1 Inhibitor
Cat. No. : RP19989
Lot. No. : PE7041710
Theoretical MW : 459.59
Observed MW : 459.40

Interface : ESI
Nebulizing Gas Flow : 1.5L/min
CDL Temp : 250°C
Block Temp : 200°C

Interface Bias
Drying Gas Flow
T.Flow
B.conc

: +4.5 kV
: 5 L/min
: 0.2 ml/min
: 50%H₂O/50%MeOH