


Mass Spectrum



Sample Information

Acquired by : Gary
 Month-Day Processed : 11/01/18
 Time Processed : 8:50:13
 Injection Volume : 0.2
 Sample Name : HIV-1 TAT (48-60)
 Cat. No. : RP20423
 Lot. No. : PE4731810
 Theoretical MW : 1719.00
 Observed MW : 1719.0

Interface : ESI
 Nebulizing Gas Flow : 1.5L/min
 CDL Temp : 250°C
 Block Temp : 200°C

Equipment : GK11010007
 Interface Bias : +4.5 kV
 Drying Gas Flow : 5 L/min
 T.Flow : 0.2 ml/min
 B.conc : 50% H₂O/50% MeOH