


Mass Spectrum



Sample Information

User : Gary
Month-Day Processed : 04/11/18
Time Processed : 10:32:37
Injection Volume : 0.2
Sample Name : β -Amyloid (1-40), rat
Cat No. : RP10016-0.5
Lot No. : PE4051804
Theoretical MW : 4233.74
Observed MW : 4234.40

Interface

Interface : ESI
Nebulizing Gas Flow : 1.5L/min
CDL Temp : 250°C
Block Temp : 200°C

Interface Bias

Interface Bias : +4.5kV
Drying Gas Flow : 5 L/min
T.Flow : 0.2ml/min
B.conc : 50%H₂O/50%MeOH