


Mass Spectrum



Sample Information

User : Gary
 Month-Day Processed : 03/04/18
 Time Processed : 14:10:55
 Injection Volume : 0.2
 Sample Name : PACAP(1-27)₃human,ovine,rat
 Cat No. : RP10334
 Lot No. : PE6051803
 Theoretical MW : 3147.62
 Observed MW : 3148.40

Interface : ESI
 Nebulizing Gas Flow : 1.5L/min
 CDL Temp : 250 °C
 Block Temp : 200 °C

Interface Bias : +4.5kV
 Drying Gas Flow : 5 L/min
 T.Flow : 0.2ml/min
 B.conc : 50%H₂O/50%MeOH