


# Mass Spectrum



Sample Information  
Month-Day Processed : 08/05/21  
Time Processed : 21:15:58  
Injection Volume : 0.2  
Sample Name : MDYKDHGDYKDHIDYKDDDDDK peptide  
Cat. No. : RP21087  
Lot. No. : E4712021G  
Theoretical MW : 2861.87  
Observed MW : 2861.6

Interface : ESI  
Nebulizing Gas Flow : 1.5L/min  
CDL Temp : 250  
Block Temp : 200

Equipment : ZJ21010035  
Interface Bias : +4.5 kV  
Drying Gas Flow : 5 L/min  
T.Flow : 0.2 ml/min  
B.conc : 50%H2O/50%MeOH